



biolas[®]

www.biolas.net

NAVIGATOR CATHETER



NAVIGATOR CATHETER

Navigator Catheter is an optical camera routable guiding catheter system that is developed for lumbar epidural endoscopy. Besides than diagnostic epiduroscopy procedure, navigator catheter's design also allows neural decompression, treatment of epidural adhesions, injection of anti-inflammatory / steroid drugs into epidural space.

Advantages


- Light and ergonomical design
- Soft a-traumatic tip
- Easy navigability
- Adjustable tip angle (90°-180°)
- Radiopaque shaft
- Double Lumen
- Wide lumens
- Low infection risk

Technical Information

Length: 30 cm

Outer Diameter: 9F/8F

Tip Angle: (90°-180°)



biolas[®]
HEALTH INC.

Fatih Sultan Mah. 2365 St. No:19
06790 - ANKARA / TURKEY - Tel: +90 312 472 1125 - +90 312 472 1126

www.biolas.net // info@biolas.net