


biolas®

REYA
HYBRID THROMBUS MANAGEMENT



REYA VENOCAT

THROMBOLYSIS + PROTECTION

www.biolas.net

REYA® VENOCAT offers a new treatment standard for reducing risk of thromboembolism during thrombolysis and thrombectomy procedure. Venocat is a complimentary product with safety and efficiency without disadvantages of current treatment options.

Advantages

1

Reduces mortality and morbidity by eliminating recurrent pulmonary embolism

2

Can be used easily with patient groups such as where anticoagulant therapy absolutely contraindicated and also patients with severe complications such as thrombocytopenia or major bleeding.

3

Minimizes complications of permanent IVC filters (fracture, penetration, migration, rupture)

4

100 % retrievable and easy to use and extremely safe

5

Can also be used as central venous catheter by it's ports and lumens

Product Specifications

Filter Diameter	Catheter Diameter	Catheter Length	Max. Duration time in body
30mm	7F	30 cm	12Days

biolas[®]
HEALTH INC.